

Bubble Column Void Fraction and Bubble Size Measurement

Alexandra E. Sigel^{1,2}, Manjil Ray^{1,2}, Rodney O. Fox^{1,3}, Alberto Passalacqua^{1,2,3}, and Theodore J. Heindel^{1,2,3}

¹Center for Multiphase Flow Research and Education

²Department of Mechanical Engineering

³Department of Chemical and Biological Engineering

Iowa State University, Ames, IA, 50011, USA



Overview


- Motivation
- Experimental setup and details
- Overall void fraction results
- Bubble size measurements
- Future work
- Conclusions

Motivation

- Bubble column bioreactors are often used because of their desirable characteristics
- Scaleup of these bioreactors continues to be a challenge faced in chemical and fermentation industries
- Part of a larger project to validate computational fluid dynamic (CFD) tools to better understand the relationships between bubble column hydrodynamics and fermentation yield



Experimental Setup

- 11 pressure transducers used to simultaneously measure pressure at various locations
- Gas is sparged through 8-arm spider sparger
- Column inner diameter of 32.1cm
- Fluid (water) static height is 10 column diameters or 3.21 meters





Sparger Details

- 8-arm spider sparger
- 26 1-mm holes in each arm arranged in 2 rows of 13 holes, with a 30° offset between rows
- Holes facing downward offset $\pm 15^\circ$ of bottom dead center
- Total open area ratio: 0.20%
- Height from bubble column base (h) is set at $h = 31$ mm but is adjustable ($31 < h < 254$ mm)



Sparger Section View




Overall Void Fraction

- Void fraction (also called volumetric gas fraction or gas holdup) determined by:


$$\varepsilon = 1 - \frac{\Delta P}{\Delta P_0}$$

- ε = void fraction between any two pressure transducers
- ΔP = average pressure drop between two pressure transducers when there is gas flow
- ΔP_0 = average pressure drop between two pressure transducers without gas flow

Overall Void Fraction Comparison



Bubble Column Void Fraction Data




Source:

Source: Su, X., and Heindel, T.J., "Modeling Gas Holdup in Gas-Liquid-Fiber Semibatch Bubble Columns," *Industrial and Engineering Chemistry Research*, 44(24) 9355-9363 (2005).


Overall Void Fraction (cont.)

- Operating in pure heterogeneous regime due to sparger geometry and bubble coalescence
- Hysteresis is observed with “fresh” water but not after the column is operated for several hours
 - “fresh” water for trials 1 and 3




Bubble Size Image Acquisition


- 30.5 cm borescope inserted into the bubble column with exterior light sources paired with high-speed camera used for image collection
 - Others have used borescopes in homogeneous flows
- Images collected at multiple axial and radial locations paired with pre-selected superficial gas velocities
 - Superficial gas velocities for collection were determined from results and industry collaborator preferences




Borescope Image Analysis Challenges




Original Video




Original Still Frame



Binarized Frame



Processed Frame




Verified Frame

- Image information: $U_g = 2$ cm/s, 60 fps
- In-house image analysis code to determine bubble sizes by detecting bubble edges
- Bubbles that are isolated are easily identified. If there are bubbles 'stacked' the code cannot currently identify the edges effectively
 - Solved with lighting to create a broader range of intensities with the images
 - Currently improving process

Next Steps

- Improve image acquisition for bubble size distribution measures
- Use our X-ray flow visualization facility to complete X-ray Computed Tomography (XCT) imaging
 - Provides time-average local void fraction within the imaging region
- Introduce a fermentation broth simulant as the working fluid to conduct overall and local void fraction and bubble size measurements



Conclusions

- Bubble column bioreactors are used in the fermentation industry but are a challenge for scale-up
- Overall void fraction data provides information across a wide range of superficial gas velocities
- Operating in the pure heterogeneous regime provides challenges for bubble size distribution data collection
- Bubbles are captured with a borescope
 - Addressing imaging challenges

Acknowledgments

- This material is based on research sponsored by the Air Force under agreement number FA8650-21-2-5028. The U.S. Government is authorized to reproduce and distribute reprints for Governmental purposes notwithstanding any copyright notation therein. This project was made possible with the support of The Bioindustrial Manufacturing and Design Ecosystem (BioMADE); the content expressed herein is that of the authors and does not necessarily reflect the views of BioMADE.
- The input from Cargill and Geno are greatly appreciated.

Bubble Column Void Fraction and Bubble Size Measurement

Alexandra E. Sigel^{1,2}, Manjil Ray^{1,2}, Rodney O. Fox^{1,3}, Alberto Passalacqua^{1,2,3}, and Theodore J. Heindel^{1,2,3}

¹Center for Multiphase Flow Research and Education

²Department of Mechanical Engineering

³Department of Chemical and Biological Engineering

Iowa State University, Ames, IA, 50011, USA

